



"We have been using the Q-Box AQUA respirometry package for a year and have been impressed with the quality and versatility of the system. The technicians at Qubit are outstanding; they were able to modify the system to suit our particular needs and they have always been quick to respond when technical support was required."

*Dr. Tyson MacCormack
Mount Alison University, New Brunswick, Canada*

Q-BOX SERIES Q-BOX AQUA AQUATIC RESPIROMETRY PACKAGE



The Q-Box AQUA is used to measure metabolic rate in aquatic animals using intermittent flow respirometry. O₂ consumption by the aquatic organism is measured with a dissolved O₂ electrode and software. Software provides calculations of the oxygen consumption rate (VO₂). Measurements are made sequentially without removing the animal from the sample chamber, and the chamber is flushed between measurements with fresh water under ambient conditions. This intermittent flow method overcomes the disadvantage of continuous closed respirometry in which the dissolved O₂ level in the chamber may decrease to a hypoxic level.

Features:


- Optical DO Probe and Water Bath
- Small Respirometry Chamber (140ml)
- Large Respirometry Chamber (660ml)
- Temperature & Salinity Probe for Water Bath
- Liquid Pump (1 LPM) & 3-Way Valve
- Digital Control Unit for 3-Way Valve Control
- 4 Channel Data Interface
- Data Acquisition Software
- Absolute Pressure Sensor
- Rugged Weatherproof Case
- Optional Battery Pack

Applications:


- Aquatic Respirometry Studies
- Vertebrate and invertebrate aquatic studies
- Lab and Field Studies



1573 John Counter Blvd., Kingston, Ontario, Canada K7M3L5
www.qubitsystems.com / info@qubitsystems.com
World Wide 01.613.384.1977
North America 888.262.2219



SOFTWARE



Sample Data from a Goldfish




During an experiment, the animal is placed in the chamber which is immersed in a water bath. Water from the bath under ambient conditions is flushed through the chamber to establish baseline measurements. Respiration measurements are made by recirculating the water through the chamber. In this recirculation phase, dissolved O₂ (DO) declines in the closed loop as the animal consumes O₂. From the rate of decline (i.e. the slope of the DO decrease), the respiratory O₂ consumption (VO₂) is calculated in the software. During the flush period, the 3-way valve is activated by the software to open the loop. The pump draws water from the bath to flush the chamber with fresh water. During the recirculation phase the 3-way valve is activated by the software to close the loop. The flush-recirculation cycle is then repeated until the experiment is terminated. The times of the flush and recirculation periods are selected by the user. For respirometry measurements of very small aquatic animals consider our Mini-AQUA Package (miniaturized version of the Q-Box AQUA).

Q-Box sensors and analyzers are compatible with any data acquisition system you may already possess.

SOFTWARE FEATURES

- 4 Data Channels
- Variable Data Acquisition Rate
- Graph, Table and Meter Displays
- Suite of Analysis Functions: (Statistics, Curve Fit, Integration etc.)
- Control of Devices using Data Input.
- Multi-Page Reporting
- Data Replay Function
- Simple Equation Input and Editing
- Numerous Mathematical Functions
- Automated Calculations On-the-Fly



Mini Aqua Chamber



Liquid Pump



1573 John Counter Blvd., Kingston, Ontario, Canada K7M3L5
 www.qubitsystems.com / info@qubitsystems.com
 World Wide 01.613.384.1977
 North America 888.262.2219

